suiteHEART[®] Software cMRI Analysis Software Release Notes suiteHEART[®]-Software cMRI-Analyse-Software Versionshinweise Logiciel suiteHEART[®] Logiciel d'analyse cMRI Notes de publication Software suiteHEART[®] Software per analisi cMRI Note di rilascio Λογισμικό suiteHEART[®] cMRI, λογισμικό ανάλυσης Σημειώσεις έκδοσης "suiteHEART[®]" programinė įranga cMRI analizės programinė įranga Išleidimo pastabos Software suiteHEART[®] Software de análisis cMRI Notas de publicación suiteHEART[®] cMRI Analysprogramvara Versionsinformation suiteHEART[®] Software cMRI Analiz Yazılımı Sürüm Notları Software-ul suiteHEART[®] Software de analiză cMRI Note privind versiunea suiteHEART[®] Software cMRI-analysesoftware Opmerkingen bij deze uitgave suiteHEART[®] Software de análise cMRI Notas de lançamento suiteHEART[®] szoftver cMRI elemzőszoftver Kibocsátási megjegyzések

NeoSoft, LLC

NS-03-041-0003 Rev. 2 Copyright 2021 NeoSoft, LLC All rights reserved

suiteHEART[®] Software cMRI Analysis Software

Release Notes

NeoSoft, LLC

NS-03-041-0003 Rev. 2 Copyright 2021 NeoSoft, LLC All rights reserved

Sales: <u>orders@neosoftmedical.com</u> Service: <u>service@neosoftmedical.com</u>

The Instructions for Use for this device are provided electronically in Portable Document Format, (.pdf). A pdf viewer is required to view the Instructions for Use. A paper copy of the Instructions for Use can be provided upon request, at no cost, within 7 calendar days, by emailing service@neosoftmedical.com.

Instructions for Use may be accessed in the following ways:

- 1. After launching the application, click "Help" or "About" from the main screen. Select the "Instructions for Use" option. The Instructions for Use will open in a pdf viewer.
- If the original installation package received from NeoSoft is available, open the zip file and navigate to the "Documentation" folder, then the "Instructions for Use" folder and double click on the Instructions for Use.pdf in your language, denoted by EN English, FR French, DE German, EL Greek, IT Italian, LT Lithuanian, ES Spanish, SV Swedish, TR Turkish, RO Romanian, NL Dutch, ZH_CN Simplified Chinese, PT-PT Portuguese, HU Hungarian.
- Navigate to the folder where the application is installed. Locate the "Documentation" folder, open the "Instructions for Use" folder and double click on the Instructions for Use.pdf in your language, denoted by EN - English, FR - French, DE - German, EL - Greek, IT - Italian, LT - Lithuanian, ES -Spanish, SV - Swedish, TR - Turkish, RO - Romanian, NL – Dutch, ZH_CN - Simplified Chinese, PT-PT - Portuguese, HU - Hungarian.
- Electronic copies of Instructions for Use are also available at <u>www.neosoftllc.com/neosoft/product_manuals/</u> for at least 15 years from last date of manufacture.

Medical device directive

This product conforms with the requirements of council directive 93/42/EEC concerning medical devices when it bears the following CE mark of Conformity:

European Representatives:

EMERGO EUROPE Prinsessegracht 20 2514 AP The Hague The Netherlands

Canada:

Health Canada device license number: 99490

Malaysian Medical Device Registration Number: GB10979720-50888

Malaysia Authorized Representative: Emergo Malaysia Sdn. Bhd. Level 16, 1 Sentral Jalan Stesen Sentral 5 KL Sentral, 50470 Kuala Lumpur MALAYSIA

CAUTION: United States federal law restricts this device to sale by, or on the order of, a physician.

Software Revision History

suiteHEART[®] Software Revision History

Software Version	Release Date	Description
1.1.1	March 2013	Initial Release
1.1.2	February 2014	Implemented software licensing options
1.1.3	April 2014	Investigational Device Release for Canada
1.1.4	May 2014	Fixed Export Cine files not retaining image manipulation settings. Fixed images not present after performing a DICOM import Corrected infarct overlay paint issue with ME Analysis Adjusted cross reference criteria Corrected slice order issue for DICOM created series Corrected Custom Series not present after reloading of an approved exam Corrected issue with suiteHEART not launching from a non-admin user
200	Ostah as 2014	account
2.0.0	October 2014	Added support for Enhanced MR Image Storage Added Valve Plane Analysis calculations Added support for Microsoft Windows 8.1 Redesigned User Interface Look and Feel Dropped support for Microsoft Windows XP
2.1.0	January 2015	Minor updates and defect fixes Added Medical Device Directive information
3.0.0	October 2015	Added Edema and Salvage mass results Added Dyssynchrony Added T1 Mapping (Research Only) Minor undates and defect fixes
3.0.1	March 2016	Added Multiuser capabilities Minor updates and defect fixes
4.0.0	October 2016	Added 3D/4D Viewer Added Nudging Tool Added Compare Mode Added Auto Phase Error Correction
4.0.1	January 2017	Minor updates and defect fixes
4.0.2	May 2017	Added tablet mode. Removed support for Microsoft Windows 8.1 Improved Auto Segmentation workflow.
4.0.3	July 2017	Added ED/ES Only mode to Auto Segmentation. Improved Auto Segmentation algorithms. Minor updates and defect fixes.
4.0.4	November 2017	Added the Viewer. Updated T1 Mapping. Added T2 Mapping. Added the Siemens Auto Series creator.
4.0.6	May 2018	Upgrade CUDA to version 9.1 Add Strain Analysis (Research Only) Add DENSE Strain Analysis (Research Only) Combine Manual/Auto Analysis modes Upgrade to JAVA 9

Software Version	Release Date	Description
4.0.7	November 2018	3D/4D Viewer Improvements Improved editing of contours Improved papillary segmentation DENSE performance improvements Automate basal line Added Undo functionality
5.0.0	July 2019	Contains updates relating to: Preprocessing, Virtual Fellow™, Series Auto Combine, Function Analysis, T2 Mapping, and general information.
5.0.1	February 2020	Contains updates relating to: Viewer, Virtual Fellow™, Function Analysis, Flow Analysis, Time Course, 3D/4D Viewer and Preferences.
5.0.2	March 2021	Added enhancements to T1 and T2 Mapping. ROI editing tool enhancements. Added additional Flow vessel categories. Enhancements to system preferences. Minor updates and defect fixes. This release is the initial software release for China.

Software Updates

T1 Mapping Analysis

It is recommended to use the scanner generated T1 map images if available for analysis.

All slice locations used for analysis should be in one series. Set the Auto Compose Series for Analysis according to vendor type in preferences Figure 1.

Figure 1: Auto Compose Series for Analysis

Siemens MyoMaps

To obtain T1 results when analyzing the original DICOM images that are similar to the generated scanner T1 map, confirm the efficiency inversion pulse used for MyoMaps MOLLI protocols. If indicated as "Non-sel IR T1 Map" on the scanner under the Contrast/Common card under Magn Preparation, the recommended inversion correction factor ICF=1.0365. For further clarification it is recommended to contact your Siemens Applications Support Specialist.

Reference: Kellman, P., Hansen, M.S. T1-mapping in the heart: accuracy and precision. *J Cardiovasc Magn Reson* 16, 2 (2014). <u>https://doi.org/10.1186/1532-429X-16-2</u>

If analyzing the original MOLLI DICOM images enter the appropriate ICF in preferences as shown in Figure 2.

- 1. Select Tools > Preferences > Edit.
- 2. Select the T1/T2 Mapping tab.
- 3. Enter the ICF according to vendor type

Figure 2: T1 Mapping Preferences

Global	Template	Macro	Print	Virtual Fel	low®	T1/T2 Mappir	ng
	Sequenc	e l	О мо	u (Satu	ration Recov	ery
	DICOM O	verlay	Non	ie 🤇	T1	ECV	
	Map Over	1ay I	O Non	ie 🤇	D T1	ECV	
		ICF					
	GE		1.0	000			
	Philip		1.0	000			
	Siem	ens	10	000			
			1.0				
	Na	ative					

Note: Users who have established normal ranges based upon the original DICOM data will need to multiply established T1 values by the recommended ICF to adjust the normal range.

T2 Mapping Analysis

It is recommended to use the scanner generated T2 map images if available for analysis.

All slice locations used for analysis should be in one series. Set the Auto Compose Series for Analysis according to vendor type in preferences Figure 1.

Selections for fitting options now include:

- Nonlinear (default)
- Linear (non-weighted)

The Nonlinear fitting algorithm no longer estimates the background noise.

Note: To obtain T2 results using Siemens original DICOM images that are similar to the generated Siemens scanner T2 map select Linear Fitting.

Viewer

If the series name is truncated in the series list a tool-tip will now appear with the series name when hovered over.

Function Analysis

Result Table

Results can be reordered by clicking and dragging the result to a new position on the table. It is recommended to define the order in preferences.

ROI Editing Tools

ROI Point Spline

Additional functionality now includes:

- Alt+ Left mouse generates a corner point.
- Clicking the first point closes the contour.
- Clicking on the contour directly generates a point.
- Delete + cursor on point removes a point.
- Dragging a point close to a neighboring point removes the neighboring point.
- If the number of points becomes less than 3, the ROI will be deleted.

Note: The point spline functionality applies to all analysis except for 3D/4D Flow Viewer. The following copy/paste and translate functionality is available in all analysis modes except for PFO analysis.

- Ctrl +C Copy ROI
- Ctrl+V Paste ROI
- Ctrl +S Smooth ROI
- Left and Right Arrow Keys W & S Keys
- Up and Down Arrow Keys A & D Keys

Flow Analysis

- When selecting regurgitant mode for a current analyzed flow series or a previous analyzed series the negative flow value will be presented as a absolute value.
- Additional editing tools are now available on the interface when the ROI is selected in edit mode. (Refer to the suiteHEART[®] Software Instructions for Use)
- The y axis for Pressure Half Time mode is now displayed as velocity (mm/s).

Release Notes

Preferences

In the multi-user environment, changes to preferences should only be performed by one user.

Templates

• The results for cardiac output (CO) for the LV and RV are reported in I/min. When entering parameters for the gender specific equation for Z-scores, the **a** parameter for both male and female must be converted to I/min. When entering parameters from reference papers the units used must be the same.

Flow Analysis

- The preference selection to show advanced options of Positive, Negative, Peak Envelope and Peak Absolute modes have been removed.
- When reviewing flow analysis from prior studies confirm the vessel category for each flow result.

Preprocessing

• Auto Detect is not supported for Canon Medical studies.

T1 Mapping

- Editing a blood pool ROI using the original DICOM images and dragging the ROI to the edge of the Endo ROI will remove the blood pool ROI from all times points. The blood pool T1 value will still be present on the table. Select the trashcan and select Blood Pool to reset.
- Legacy studies with previous T1 analysis, the analysis will need to be redone.

T2 Mapping

• Legacy studies with previous T2 analysis, the analysis will need to be redone.

Myocardial Evaluation

• For the Early Enhancement analysis if a slice classification is set to Undefined the results from that slice will still be included in the results. To remove the slice from the analysis, delete all ROIs.

3D/4D Flow Viewer

• For 4D Flow the cine option will be enabled after image loading has been have completed.

General

- When using suiteHEART[®] Software at a screen resolution higher than 1920x1080 the matrix remains at the same size as it appears at 1920x1080.
- When networking a study to suiteDXT, do not perform any analysis until the study has been transferred completely.

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