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Sales: [orders@neosoftmedical.com](mailto:orders@neosoftmedical.com)

Service: [service@neosoftmedical.com](mailto:service@neosoftmedical.com)

The Instructions for Use for this device are provided electronically in Portable Document Format, (.pdf). A pdf viewer is required to view the Instructions for Use. A paper copy of the Instructions for Use can be provided upon request, at no cost, within 7 calendar days, by emailing [service@neosoftmedical.com](mailto:service@neosoftmedical.com).

Instructions for Use may be accessed in the following ways:

1. After launching the application, click “Help” or “About” from the main screen. Select the “Instructions for Use” option. The Instructions for Use will open in a pdf viewer.
2. If the original installation package received from NeoSoft is available, open the zip file and navigate to the “Documentation” folder, then the “Instructions for Use” folder and double click on the Instructions for Use.pdf in your language, denoted by EN - English, FR - French, DE - German, EL - Greek, IT - Italian, LT - Lithuanian, ES - Spanish, SV - Swedish, TR - Turkish, RO - Romanian, NL - Dutch, ZH\_CN - Simplified Chinese, PT-PT - Portuguese, HU - Hungarian, JA - Japanese, VI - Vietnamese.
3. Navigate to the folder where the application is installed. Locate the “Documentation” folder, open the “Instructions for Use” folder and double click on the Instructions for Use.pdf in your language, denoted by EN - English, FR - French, DE - German, EL - Greek, IT - Italian, LT - Lithuanian, ES - Spanish, SV - Swe-dish, TR - Turkish, RO - Romanian, NL – Dutch, ZH\_CN - Simplified Chinese, PT-PT - Portuguese, HU - Hungarian, JA - Japanese, VI - Vietnamese.
4. Electronic copies of Instructions for Use are also available at [www.neosoftllc.com/neosoft/product\\_manuals/](http://www.neosoftllc.com/neosoft/product_manuals/) for at least 15 years from last date of manufacture.

# Medical device directive

This product conforms with the requirements of council directive 93/42/EEC concerning medical devices when it bears the following CE mark of Conformity:

This product is a Medical Device:

## European Representatives:

EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

## EU Importer:

MedEnvoy Global B.V.  
Pr. Margrietplantsoen 33, Suite 123,  
2595 AM The Hague,  
The Netherlands

## Swiss Authorized Representative and Importer:

MedEnvoy Switzerland  
Gotthardstrasse 28  
6300 Zug  
Switzerland

## Malaysia Authorized Representative:

Emergo Malaysia Sdn. Bhd.  
Level 16, 1 Sentral Jalan Stesen Sentral 5 KL  
Sentral, 50470  
Kuala Lumpur MALAYSIA

Malaysian Medical Device Registration Number: GB10979720-50888

Canada: Health Canada device license number: 99490

**CAUTION:** *United States federal law restricts this device to sale by, or on the order of, a physician.*

Any serious incident that has occurred in relation to this device should be reported to NeoSoft and the competent authority of your Member State.

# Software Revision History

## suiteHEART® Software Revision History

| Software Version | Release Date  | Description                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------------------|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1.1.1            | March 2013    | Initial Release                                                                                                                                                                                                                                                                                                                                                                                                                              |
| 1.1.2            | February 2014 | Implemented software licensing options                                                                                                                                                                                                                                                                                                                                                                                                       |
| 1.1.3            | April 2014    | Investigational Device Release for Canada                                                                                                                                                                                                                                                                                                                                                                                                    |
| 1.1.4            | May 2014      | Fixed Export Cine files not retaining image manipulation settings. Fixed images not present after performing a DICOM import<br>Corrected infarct overlay paint issue with ME Analysis<br>Adjusted cross reference criteria<br>Corrected slice order issue for DICOM created series<br>Corrected Custom Series not present after reloading of an approved exam<br>Corrected issue with suiteHEART not launching from a non-admin user account |
| 2.0.0            | October 2014  | Translated application to German and French<br>Added support for Enhanced MR Image Storage<br>Added Valve Plane Analysis calculations<br>Added support for Microsoft Windows 8.1<br>Redesigned User Interface Look and Feel<br>Dropped support for Microsoft Windows XP                                                                                                                                                                      |
| 2.1.0            | January 2015  | Minor updates and defect fixes<br>Added Medical Device Directive information                                                                                                                                                                                                                                                                                                                                                                 |
| 3.0.0            | October 2015  | Added Edema and Salvage mass results<br>Added Dyssynchrony<br>Added T1 Mapping (Research Only)<br>Minor updates and defect fixes                                                                                                                                                                                                                                                                                                             |
| 3.0.1            | March 2016    | Added Multiuser capabilities<br>Minor updates and defect fixes                                                                                                                                                                                                                                                                                                                                                                               |
| 4.0.0            | October 2016  | Added 3D/4D Viewer<br>Added Nudging Tool<br>Added Compare Mode<br>Added Auto Phase Error Correction                                                                                                                                                                                                                                                                                                                                          |
| 4.0.1            | January 2017  | Minor updates and defect fixes                                                                                                                                                                                                                                                                                                                                                                                                               |
| 4.0.2            | May 2017      | Added tablet mode.<br>Removed support for Microsoft Windows 8.1<br>Improved Auto Segmentation workflow.                                                                                                                                                                                                                                                                                                                                      |
| 4.0.3            | July 2017     | Added ED/ES Only mode to Auto Segmentation.<br>Improved Auto Segmentation algorithms.<br>Minor updates and defect fixes.                                                                                                                                                                                                                                                                                                                     |
| 4.0.4            | November 2017 | Added the Viewer.<br>Updated T1 Mapping.<br>Added T2 Mapping.<br>Added the Siemens Auto Series creator.                                                                                                                                                                                                                                                                                                                                      |
| 4.0.6            | May 2018      | Upgrade CUDA to version 9.1<br>Add Strain Analysis (Research Only)<br>Add DENSE Strain Analysis (Research Only)<br>Combine Manual/Auto Analysis modes<br>Upgrade to JAVA 9                                                                                                                                                                                                                                                                   |
| 4.0.7            | November 2018 | 3D/4D Viewer Improvements<br>Improved editing of contours<br>Improved papillary segmentation<br>DENSE performance improvements<br>Automate basal line<br>Added Undo functionality                                                                                                                                                                                                                                                            |
| 5.0.0            | July 2019     | Contains updates relating to: Preprocessing, Virtual Fellow™, Series Auto Combine, Function Analysis, T2 Mapping, and general information.                                                                                                                                                                                                                                                                                                   |

| Software Version | Release Date   | Description                                                                                                                                                                                                                                  |
|------------------|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 5.0.1            | February 2020  | Contains updates relating to: Viewer, Virtual Fellow™, Function Analysis, Flow Analysis, Time Course, 3D/4D Viewer and Preferences.                                                                                                          |
| 5.0.2            | March 2021     | Added enhancements to T1 and T2 Mapping. ROI editing tool enhancements. Added additional Flow vessel categories. Enhancements to system preferences. Minor updates and defect fixes. This release is the initial software release for China. |
| 5.0.2            | June 2021      | Added MD symbol, EU importer reference, incident reporting information.                                                                                                                                                                      |
| 5.0.3            | September 2021 | Enhancements including: Auto Update, Preprocessing for T1 and T2 Mapping, Function editing, system preferences, exporting to XLS and XML. Minor updates and defect fixes. This is the initial software release for Japanese and Vietnamese.  |

# Software Updates

## Auto Update

The Auto Update feature allows for the launching of a study with background processing. Refer to the suiteHEART® Instruction for Use and the suiteDXT Instructions for Use.

## Viewer

The velocity color overlay is now supported for phase contrast images to provide quick identification. The opacity can be adjusted in preferences. To remove the color overlay set the opacity to 0%. **Tools->Preferences->Edit->Global->Viewer**

## Quick Keys

- Pressing the *Ctrl key + middle mouse button* allows for the scrolling through phases.
- The Region of Interest tool now reports the pixel size of the ROI.
- Quick Keys are available for generic annotations.


|                    |         |
|--------------------|---------|
| Linear             | Shift+1 |
| Crosshair          | Shift+2 |
| Region of Interest | Shift+3 |
| Label              | Shift+4 |
| Angle              | Shift+5 |

## Virtual Fellow®

Linear measurements can now be performed for the long axis viewports.

## T1 and T2 Mapping Analysis

- Preprocessing is now supported for T1 and T2 Mapping. The series type of either Time Series or Map is supported. To enable the preprocessing and set and the selection of the series type, from the suiteDXT application, **Configuration-> Application Launch tab-> T1 Mapping or T2 Mapping.**
- In suiteHEART® Software, adjust the endo and epi contours using the Endo (red) or Epi (green) offset control bar.

## Function Analysis


- For Regional Analysis the default ranges for Thickness, Pct. Thickening and Wall Motion can be set in preferences. **Tools > Preferences > Edit > Global > Function**
- Auto segmentation for the LA and RA is supported for Strain Analysis\*. Please refer to the suiteHEART Instruction for Use
- Displayed in the Editor viewport are selections for toggling between the two edit modes.

| Tool                                                                              | Description  |
|-----------------------------------------------------------------------------------|--------------|
|  | Limit ROI    |
|  | No Limit ROI |

- Triangle indicators have been added to the editor viewport and volume curve to indicate end-diastole (red) and end-systole (blue).
- Editing mode now persists for endo and epi contours. **Tools->Preferences->Edit->Global->Function** check the selection for **Persist Endo/Epi ROI edit mode**.

## Flow Analysis

- The ml/beat and l/min selections are available from the file pull-down menu as shown below.
- The selection for Positive Flow (calculated flow above the baseline) has been added to the Flow Tab interface as shown.

## 3D/4D Flow Viewer

1x1 layout. Double click on any 2x2 viewport to toggle the layout.

## Structured Reporting

Upon the approving of an exam, automatic export of results can be performed to either Acuity, Powerscribe or Precession.

## Preferences

12-Hour or 24-Hour format is supported. **Tools->Preferences->Edit->Global->General**

# Release Notes

## Operating System Support

suiteHEART® Software version 5.0.3 will be the last version of software supported on Windows 7.

## Flow Analysis

When reviewing flow analysis from prior studies confirm the vessel category for each flow result.

## Myocardial Evaluation

When performing T2 analysis for edema do not manually delete a single normal ROI. Delete all normal ROIs.

## T1 Mapping

- Editing a blood pool ROI using the original DICOM images and dragging the ROI to the edge of the Endo ROI will remove the blood pool ROI from all times points. The blood pool T1 value will still be present on the table. Select the trashcan and select Blood Pool to reset.
- When selecting the ECV color overlay, if the positioning of the anatomy in the FOV has changed from native to post, the ROIs will be misaligned. Results are not affected.

## 3D/4D Flow Viewer

- If a study has been launched and 4D flow results have been generated using ROI 1 for a vessel category, if additional 2D phase contrast series are then added and being networked and auto update performed the 2D flow series will not be preprocessed.
- For 4D Flow the cine option will be enabled after image loading has been have completed.
- If using single oblique or double oblique mode when using the locate feature the oblique lines may not represent the location of the anatomy in the viewport. Reselect the oblique tool to reset.

## Preferences

- In the multi-user environment, changes to preferences should only be performed by one user.

## Templates

- The results for cardiac output (CO) for the LV and RV are reported in l/min. When entering parameters for the gender specific equation for Z-scores, the a parameter for both male and female must be converted to l/min. When entering parameters from reference papers the units used must be the same.

## General

- When using suiteHEART at a screen resolution higher than 1920x1080 the matrix remains at the same size as it appears at 1920x1080.
- For preprocessing Auto Detect is not supported for Canon Medical studies.

## Strain Analysis\*

- It is recommended to perform ventricular strain analysis first and atrial strain analysis second.
- Due to ROI overlap there may be instances when selecting the LV Epi ROI, the strain mode will switch from LV results to RV results.
- Strain results are based upon smooth contours, if contours are added manually, create the contours using the smooth mode in function analysis.

***\*Strain Analysis is only for research purposes. Claims have not been evaluated by the FDA. Use pursuant to company instruction and research agreement. NeoSoft, LLC proprietary and confidential.***

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All rights reserved