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NEOSOFT

NeoSoft, LLC

NS-03-043-0011 Rev. 1
Copyright 2022 NeoSoft, LLC
All rights reserved

suiteHEART® Software

cMRI Analysis Software

Release Notes

NeoSoft, LLC

NS-03-043-0011-EN Rev. 1
Copyright 2022 NeoSoft, LLC
All rights reserved

Manufacturer

NeoSoft, LLC
N27 W23910A Paul Road
Pewaukee, WI 53072 USA

Phone: 262-522-6120
website: www.neosoftllc.com

Sales: orders@neosoftmedical.com
Service: service@neosoftmedical.com

The Instructions for Use for this device are provided electronically in Portable Document Format, (.pdf). A pdf viewer is required to view the Instructions for Use. A paper copy of the Instructions for Use can be provided upon request, at no cost, within 7 calendar days, by emailing service@neosoftmedical.com.

Instructions for Use may be accessed in the following ways:

1. After launching the application, click “Help” or “About” from the main screen. Select the “Instructions for Use” option. The Instructions for Use will open in a pdf viewer.
2. If the original installation package received from NeoSoft is available, open the zip file and navigate to the “Documentation” folder, then the “Instructions for Use” folder and double click on the Instructions for Use.pdf in your language, denoted by EN - English, FR - French, DE - German, EL - Greek, IT - Italian, LT - Lithuanian, ES - Spanish, SV - Swedish, TR - Turkish, RO - Romanian, NL - Dutch, PT-PT - European Portuguese, HU - Hungarian, JA - Japanese, VI - Vietnamese, ZH-CN - Simplified Chinese, PT-BR - Brazilian Portuguese.
3. Navigate to the folder where the application is installed. Locate the “Documentation” folder, open the “Instructions for Use” folder and double click on the Instructions for Use.pdf in your language, denoted by EN - English, FR - French, DE - German, EL - Greek, IT - Italian, LT - Lithuanian, ES - Spanish, SV - Swedish, TR - Turkish, RO - Romanian, NL - Dutch, PT-PT - European Portuguese, HU - Hungarian, JA - Japanese, VI - Vietnamese, ZH-CN - Simplified Chinese, PT-BR - Brazilian Portuguese.
4. Electronic copies of Instructions for Use are also available at www.neosoftllc.com/neosoft/product_manuals/ for at least 15 years from last date of manufacture.
5. suiteHEART / suiteDXT Software Group Basic UDI DI – 850001088039.

Medical device directive

This product conforms with the requirements of council directive 93/42/EEC concerning medical devices when it bears the following CE mark of Conformity:

This product is a Medical Device:	Importer:

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Malaysia Authorized Representative:

Emergo Malaysia Sdn. Bhd.
Level 16, 1 Sentral Jalan Stesen Sentral 5 KL Sen-
tral, 50470
Kuala Lumpur MALAYSIA

Malaysian Medical Device Registration Number:
GB10979720-50888

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Importer / Registration Holder in Brazil:

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Canada: Health Canada device license number: 99490

suiteHEART 5.1.0 UDI: 00850001088015510

CAUTION: *United States federal law restricts this device to sale by, or on the order of, a physician.*

Any serious incident that has occurred in relation to this device should be reported to NeoSoft and the competent authority of your Member State.

suiteHEART® Software Revision History

Software Version	Release Date	Description
1.1.1	March 2013	Initial Release
1.1.2	February 2014	Implemented software licensing options
1.1.3	April 2014	Investigational Device Release for Canada
1.1.4	May 2014	Fixed Export Cine files not retaining image manipulation settings. Fixed images not present after performing a DICOM import Corrected infarct overlay paint issue with ME Analysis Adjusted cross reference criteria Corrected slice order issue for DICOM created series Corrected Custom Series not present after reloading of an approved exam Corrected issue with suiteHEART not launching from a non-admin user account
2.0.0	October 2014	Translated application to German and French Added support for Enhanced MR Image Storage Added Valve Plane Analysis calculations Added support for Microsoft Windows 8.1 Redesigned User Interface Look and Feel Dropped support for Microsoft Windows XP
2.1.0	January 2015	Minor updates and defect fixes Added Medical Device Directive information
3.0.0	October 2015	Added Edema and Salvage mass results Added Dyssynchrony Added T1 Mapping (Research Only) Minor updates and defect fixes
3.0.1	March 2016	Added Multiuser capabilities Minor updates and defect fixes
4.0.0	October 2016	Added 3D/4D Viewer Added Nudging Tool Added Compare Mode Added Auto Phase Error Correction
4.0.1	January 2017	Minor updates and defect fixes
4.0.2	May 2017	Added tablet mode. Removed support for Microsoft Windows 8.1 Improved Auto Segmentation workflow
4.0.3	July 2017	Added ED/ES Only mode to Auto Segmentation. Improved Auto Segmentation algorithms. Minor updates and defect fixes
4.0.4	November 2017	Added the Viewer Updated T1 Mapping Added T2 Mapping Added the Siemens Auto Series creator
4.0.6	May 2018	Upgrade CUDA to version 9.1 Add Strain Analysis (Research Only) Add DENSE Strain Analysis (Research Only) Combine Manual/Auto Analysis modes Upgrade to JAVA 9

Software Version	Release Date	Description
4.0.7	November 2018	3D/4D Viewer Improvements Improved editing of contours Improved papillary segmentation DENSE performance improvements Automate basal line Added Undo functionality
5.0.0	July 2019	Contains updates relating to: Preprocessing, Virtual Fellow™, Series Auto Combine, Function Analysis, T2 Mapping, and general information
5.0.1	February 2020	Contains updates relating to: Viewer, Virtual Fellow™, Function Analysis, Flow Analysis, Time Course, 3D/4D Viewer and Preferences
5.0.2	March 2021	Added enhancements to T1 and T2 Mapping. ROI editing tool enhancements. Added additional Flow vessel categories. Enhancements to system preferences. Minor updates and defect fixes. This release is the initial software release for China.
5.0.2	June 2021	Added MD symbol, EU importer reference, incident reporting information
5.0.3	September 2021	Enhancements including: Auto Update, Preprocessing for T1 and T2 Mapping, Function editing, system preferences, exporting to XLS and XML. Minor updates and defect fixes. This is the initial software release for Japanese and Vietnamese.
5.0.4	November 2021	Minor updates and defect fixes
5.1.0	December 2022	Added support for Microsoft Windows 11 and removed support for Microsoft Windows 7. Added Menu Driven Reporting interface, as well as other feature updates. This is the initial software release for Brazilian Portuguese.

Software Updates

Function Analysis

- The interface now has three cross reference viewports with the ability to scroll.
- Custom measurements under Other can now be reordered.
- Myocardial Contraction Fraction (MCF) result has been added.
- Information about the updated matrix view functionality can be found in the suiteHEART Instructions for Use.
- Auto combine for long axis views is available and enabled as default. (Note: Preprocessing must be performed with Virtual Fellow.)

Flow Analysis

The velocity color overlay opacity can be set in Preferences as a default display.

3D/4D Flow Viewer

- Scroll capability in viewports is available by using the middle mouse button.
- The new Quick key: Shift+1 applies linear annotations.
- Use Shift +Ctrl for target cursor.

Menu Driven Reporting

The Report View can be access by clicking  in the lower right had corner of the interface or by performing Alt+R. To facilitate the reporting of cardiac images, two monitors are recommended.

Refer to the suiteHEART Instructions for Use for more information on the updates.

T1/T2 Mapping

- Endo/Epi offsets can now be performed on a single slice.
- T1/T2 normal ranges for templates are now able to be specified for 1.5T and 3T.
- Copy and paste for local ROI is now available.

T2* Analysis

T2* normal ranges for templates can now be specified for 1.5T and 3T.

Viewer

Expand all series by clicking + (the plus key) and collapse using - (the minus key).

Release Notes

Strain Analysis*

Previous analysis performed on prior software releases for atrial strain will not be valid and will need to be reanalyzed.

*Strain Analysis is only for research purposes. Claims have not been evaluated by the FDA. Use pursuant to company instruction and research agreement. NeoSoft, LLC proprietary and confidential.

Flow Analysis

When reviewing flow analysis from prior studies confirm the vessel category for each flow result.

Myocardial Evaluation

When performing T2 analysis for edema do not manually delete a single normal ROI; delete all normal ROIs.

T1 Mapping

Editing a blood pool ROI using the original DICOM images and dragging the ROI to the edge of the Endo ROI will remove the blood pool ROI from all time points. The blood pool T1 value will still be present on the table. Select the trashcan and select Blood Pool to reset.

When selecting the ECV color overlay, if the positioning of the anatomy in the FOV has changed from native to post the ROIs will be misaligned. Results are not affected.

3D/4D Flow Viewer

If a study has been launched and 4D flow results have been generated using ROI 1 for a vessel category, then if additional 2D phase contrast series are added and being networked and auto update is performed, the 2D flow series will not be preprocessed.

For 4D Flow the cine option will be enabled after image loading has been completed.

If using single or double oblique mode when using the locate feature, the oblique lines may not represent the location of the anatomy in the viewport. Reselect the oblique tool to reset.

Preferences

In the multi-user environment, changes to preferences should only be performed by one user.

Templates

The results for cardiac output (CO) for the LV and RV are reported in l/min. When entering parameters for the gender specific equation for Z-scores, then a parameter for both male and female must be converted to l/min. When entering parameters from reference papers the units used must be the same.

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All rights reserved (無断複写・複製・転載禁止)