suiteHEART® Software cMRI Analysis Software Release Notes
suiteHEART®-Software cMRI-Analyse-Software Versionshinweise
Logiciel suiteHEART® Logiciel d'analyse cMRI Notes de publication
Software suiteHEART® Software per analisi cMRI Note di rilascio
Λογισμικό suiteHEART® cMRI, λογισμικό ανάλυσης Σημειώσεις έκδοσης
„suiteHEART® “ programinė įranga cMRI analizės programinė įranga Leidimo pastabos
Software suiteHEART® Software de análisis cMRI Notas de publicación
suiteHEART® cMRI Analysprogramvara Versionsinformation
suiteHEART® Yazılımı cMRI Analiz Yazılımı Sürüm Notları
Software-ul suiteHEART® Software de analiză cMRI Note privind versiunea
suiteHEART® Software cMRI-analysesoftware Opmerkingen bij deze uitgave
Software suiteHEART® Software de análise cMRI Notas de lançamento
suiteHEART® szoftver cMRI elemzőszoftver Kibocsátási megjegyzések
Phần mềm suiteHEART® Phần mềm phân tích cMRI Thông tin phiên bản
Software suiteHEART® Software de análise cMRI Notas de versões
suiteHEART® ソフトウェア cMRI 解析ソフトウェア リリースノート
suiteHEART® 软件 cMRI 分析软件 发布说明
Tarkvara suiteHEART® cMRI analüüsi tarkvara Väljalaskemärkmed

NEOSOFT

NeoSoft, LLC

NS-03-046-0006 Rev. 2
Copyright 2025 NeoSoft, LLC
All rights reserved

suiteHEART® Software

cMRI Analysis Software

Release Notes

NeoSoft, LLC

NS-03-046-0006-EN Rev. 2
Copyright 2025 NeoSoft, LLC
All rights reserved

Manufacturer

NeoSoft, LLC
N27 W23910A Paul Road
Pewaukee, WI 53072 USA

Phone: 262-522-6120
website: www.neosoftllc.com

Sales: orders@neosoftmedical.com
Service: service@neosoftmedical.com

To view compliance information (Authorized Representative, Importer, Registration information) after launching the application, click “Help” or “About” from the main screen. Select the “Regulatory Information” option. The document will open in a pdf viewer.

The Instructions for Use for this device are provided electronically in Portable Document Format, (.pdf). A pdf viewer is required to view the Instructions for Use. A paper copy of the Instructions for Use can be provided upon request, at no cost, within 7 calendar days, by emailing service@neosoftmedical.com.

Instructions for Use may be accessed in the following ways:

1. After launching the application, click “Help” or “About” from the main screen. Select the “Instructions for Use” option. The Instructions for Use will open in a pdf viewer.
2. If the original installation package received from NeoSoft is available, open the zip file and navigate to the “Documentation” folder, then the “Instructions for Use” folder and double click on the Instructions for Use.pdf in your language, denoted by EN - English, FR - French, DE - German, EL - Greek, IT - Italian, LT - Lithuanian, ES - Spanish, SV - Swedish, TR - Turkish, RO - Romanian, NL - Dutch, PT-PT - European Portuguese, HU - Hungarian, JA - Japanese, VI - Vietnamese, ZH-CN - Simplified Chinese, PT-BR - Brazilian Portuguese, ET - Estonian.
3. Navigate to the folder where the application is installed. Locate the “Documentation” folder, open the “Instructions for Use” folder and double click on the Instructions for Use.pdf in your language, denoted by EN - English, FR - French, DE - German, EL - Greek, IT - Italian, LT - Lithuanian, ES - Spanish, SV - Swedish, TR - Turkish, RO - Romanian, NL - Dutch, PT-PT - European Portuguese, HU - Hungarian, JA - Japanese, VI - Vietnamese, ZH-CN - Simplified Chinese, PT-BR - Brazilian Portuguese, ET - Estonian.
4. Electronic copies of Instructions for Use are also available at www.neosoftllc.com/neosoft/product_manuals/ for at least 15 years from last date of manufacture.
5. suiteHEART / suiteDXT Software Group Basic UDI DI – 00850001088039.

suiteHEART® Software Revision History

Software Version	Release Date	Description
1.1.1	March 2013	Initial Release
1.1.2	February 2014	Implemented software licensing options
1.1.3	April 2014	Investigational Device Release for Canada
1.1.4	May 2014	Fixed Export Cine files not retaining image manipulation settings. Fixed images not present after performing a DICOM import Corrected infarct overlay paint issue with ME Analysis Adjusted cross reference criteria Corrected slice order issue for DICOM created series Corrected Custom Series not present after reloading of an approved exam Corrected issue with suiteHEART not launching from a non-admin user account
2.0.0	October 2014	Translated application to German and French Added support for Enhanced MR Image Storage Added Valve Plane Analysis calculations Added support for Microsoft Windows 8.1 Redesigned User Interface Look and Feel Dropped support for Microsoft Windows XP
2.1.0	January 2015	Minor updates and defect fixes Added Medical Device Directive information
3.0.0	October 2015	Added Edema and Salvage mass results Added Dyssynchrony Added T1 Mapping (Research Only) Minor updates and defect fixes
3.0.1	March 2016	Added Multiuser capabilities Minor updates and defect fixes
4.0.0	October 2016	Added 3D/4D Viewer Added Nudging Tool Added Compare Mode Added Auto Phase Error Correction
4.0.1	January 2017	Minor updates and defect fixes
4.0.2	May 2017	Added tablet mode. Removed support for Microsoft Windows 8.1 Improved Auto Segmentation workflow
4.0.3	July 2017	Added ED/ES Only mode to Auto Segmentation. Improved Auto Segmentation algorithms. Minor updates and defect fixes
4.0.4	November 2017	Added the Viewer Updated T1 Mapping Added T2 Mapping Added the Siemens Auto Series creator
4.0.6	May 2018	Upgrade CUDA to version 9.1 Add Strain Analysis (Research Only) Add DENSE Strain Analysis (Research Only) Combine Manual/Auto Analysis modes Upgrade to JAVA 9

Software Version	Release Date	Description
4.0.7	November 2018	3D/4D Viewer Improvements Improved editing of contours Improved papillary segmentation DENSE performance improvements Automate basal line Added Undo functionality
5.0.0	July 2019	Contains updates relating to: Preprocessing, Virtual Fellow™, Series Auto Combine, Function Analysis, T2 Mapping, and general information
5.0.1	February 2020	Contains updates relating to: Viewer, Virtual Fellow™, Function Analysis, Flow Analysis, Time Course, 3D/4D Viewer and Preferences
5.0.2	March 2021	Added enhancements to T1 and T2 Mapping. ROI editing tool enhancements. Added additional Flow vessel categories. Enhancements to system preferences. Minor updates and defect fixes. This release is the initial software release for China.
5.0.2	June 2021	Added MD symbol, EU importer reference, incident reporting information
5.0.3	September 2021	Enhancements including: Auto Update, Preprocessing for T1 and T2 Mapping, Function editing, system preferences, exporting to XLS and XML. Minor updates and defect fixes. This is the initial software release for Japanese and Vietnamese.
5.0.4	November 2021	Minor updates and defect fixes
5.1.0	December 2022	Added support for Microsoft Windows 11 and removed support for Microsoft Windows 7. Added Menu Driven Reporting interface, as well as other feature updates. This is the initial software release for Brazilian Portuguese.
5.1.1	October 2023	Added interfaces to Epic and CMR Coop. Updated 4D Flow analysis functionality. Added Viewer Image/series filter functionality.suiteHEART version 5.1.1 will be the last version of software supported on Server 2012.
5.1.2	June 2025	Increased functionality and improved enhancements.
5.1.3	November 2025	Modified Global T1 calculation and reinstated Cardiac Output measurement. Added statement regarding end of support for Windows 10 and Server 2016. Minor defect fixes.

Software Updates

Refer to the suiteHEART Instructions for Use for additional details on these features.

suiteHEART 5.1.3 Release

This is the last release that will support Microsoft® Windows 10 and Windows 2016 Server.

Function

The CO Index value is now available.

T1/T2 Mapping

The issue described in the Advisory Notice that was issued on 18Sep2025 has been resolved. To decrease the artifact of shading due to variations of signal intensity from surface coil and to improve the accuracy, the Global T1 result for Time series is now calculated segment by segment and averaged.

Reporting

Additional save report format outputs are now supported.

suiteHEART 5.1.2 Release

General

The application allows individual users to configure a subset of preferences.

Threshold editing tool for LV endocardium now supported. Select  then press and hold the Alt Key.

To format the interface into multiple displays, at the top of the interface select  .

Reporting

Approval of reports can only be performed in the report preview window.

Templates

Z-scores are now available as macros.

Viewer

Split screen compare mode is now supported by selecting .

Use the preference setting for disabling cine and series expansion.

There is a new export composer for the creation of cine/image file types for images, graphs and polar plots and DICOM files.

Define images as favorites by right mouse clicking on the viewport and selecting .

Filter for series with analysis by selecting .

Function Analysis

Automatic placement of the max wall thickness measurement.

Support for multi heartbeat real-time short axis function acquisitions.

Functional analysis is now supported for multiple series. The results that are present on the report reflect the current series selected under functional analysis.

Calculation of MAPSE and TAPSE.

Calculation of e' from the 4-chamber long axis is available. For previously analyzed studies place the MV annulus line manually to obtain e' .

The E/A curve can be displayed.

For the CAS-SCPT_SERIES_AUTO_FUNCTION_LAX series the 3-ch view will no longer be used for the LV volume calculation, only the 2- and 4-chamber will be used. Method will be Biplane Simpsons rule using only the 2Ch and 4Ch segmentation.

Atrial volumetrics have been expanded, therefore any previous macros using atrial fast results will need to be updated.

Fractional Area Change (FAC) calculation from the 4-chamber long axis view for the RV, the RV volume results have been removed.

Calculation of the CMR PCWP is only applicable in certain types of heart failure, (refer to the recommended reference for more details).

ED and ES points on LV or RV volume curves can be enlarged by hovering the cursor over the points.

Basal interpolation can be turned “On” or “Off” by selecting  .

Flow

Automatic detection and correction for aliasing both 2D and 4D.

RF Denominator can be edited without LVSV or RVSV result present.

ROIs will report only area and number of pixels.

Myocardial Evaluation

Auto detection of myocardial segments with increased signal intensity and probability thresholding is available.

Admin can change the default analysis method FWHM setting in preferences.

T1/T2 Mapping

The result table has been improved.

Motion correction for T1 Mapping time series.

Support is available for GE FIESTA and FSPGR T2 Mapping time sequences.

Slice classification color overlay is now enabled.

T2*

Auto segmentation of the ventricular septum.

Myocardial Perfusion

Preprocessing for time series.

RV insertion point is now automatically added when a LV epicardium is created.

There is a new Endo/Epi Offset editing tool.

Motion correction for time series.

3D/4D Flow Viewer

Auto segmentation and measurement reporting for 3D and auto inline flow analysis for 4D flow.

Pathlines visualization is now supported for 4D Flow.

Automatic 3D Vessel Segmentation is optimized for bSSFP sequences but supports contrast enhanced 3D MRA and contrast enhanced dual-echo water reconstructed image types.

Strain Analysis*

Refer to the Strain Analysis Instructions for Use for additional details on these features.

Strain values can be batched exported from the Data Export option on suiteDXT.

Only circumferential strain results for the RV are shown, place ED manual epicardial traces for radial strain results.

Metric values for the global myocardium can be edited.

Enhancements have been made to displayed results for short axis, long axis and atrial.

The results for Displacement and Velocity now reflect the longitudinal vector of a pixel in myocardium points to base, no matter where the pixel is located.

Quantitative Perfusion*

Refer to the Quantitative Myocardial Instructions for Use for additional details on these features.

Parametric maps generated by third party software which are grayscale or grayscale with LUT (which are converted to RGB) are supported for segmentation.

Maps with just RGB pixel data are not supported for segmentation.

***Strain and Quantitative Perfusion Analysis is only for research purposes. Claims have not been evaluated by the FDA. Use pursuant to company instruction and research agreement. NeoSoft, LLC proprietary and confidential.**

Release Notes

General

Image series with the modality set as PR are not viewable in suiteHEART, only MR is supported.

Preprocessing for 3D or 4D series for Cannon is not supported.

For the report to be saved in the report database the age must be present.

Viewer

When switching to the viewer, any series created with the Save DICOM Series option or 3D/4D created series will require that the new series be selected for viewing.

Export Composer

If multiple studies are launched the exporting of images, graphs and polar plots of different studies is allowed.

Neither the multi-slice selection nor drag and drop are supported when exporting from analysis mode or review mode.

Function Analysis

Both short axis and long axis views should have the same number of phases for segmentation.

If a template has been used with rough segmentation, the template ranges will need to be entered for rough mode. When launching previous studies, from a prior release, the template will need to be selected.

When restoring a previous study the e', MAPSE and TAPSE results will not be present if an annulus line has not been placed. For best results, place the MV and TV annulus lines on the 4-chamber cine view.

If reviewing regional wall motion analysis results while simultaneously launching more than one study, the curves will not display. Either switch analysis mode or launch a single study.

Flow Analysis

Users are responsible for placing flow ROIs on all phases across the cardiac cycle for accurate flow per beat values for 2D and 4D. If there is missing segmentation an error  will be indicated on the flow graph.

When using the aliasing auto detection and correction for 2D phase contrast, aliasing may be detected due to noise along the edge of the vessel.

When reviewing flow analysis from prior studies confirm the vessel category for each flow result.

T1/T2 Mapping

GE post processed T1 and T2 Map series should not be networked to suiteDXT, only the time series is valid for analysis.

For GE T2 Mapping it is recommended to use the Royal color overlay.

If multiple T2 Map sequence types are present in the study the preprocessing will select the most recent series for analysis.

T1 mapping will not be preprocessed if the case has been anonymized with the removal of the study date.

If the short axis slices are not exactly parallel the auto-combine will fail. Use the viewer to create a series for analysis.

When selecting the ECV color overlay, if the positioning of the anatomy in the FOV has changed from native to post, the ROIs will be misaligned. Table ECV results are not affected.

Myocardial Perfusion

If the study has only one perfusion acquisition present the default label will be Stress.

If using the MOCO for preprocessing, the original MOCO series from the scanner will not be auto combined.

3D/4D Flow Viewer

Save DICOM Series feature for Cannon 4D Flow studies is not supported.

If using single oblique or double oblique mode when using the locate feature the oblique lines may not represent the location of the anatomy in the viewport. Reselect the oblique tool to reset.

Preferences

Admin privileges are determined by your IT department. Whether you are a user with access to suiteDXT admin and suiteDXT non-admin log in depends on the install.

Macros

Macros created with fast atrial results from the prior release will need to be updated.

Templates

Normal ranges entered for templates and the categorical ranges reporting feature should match. It is not recommended to have two different range values for reporting.

The results for cardiac output (CO) for the LV and RV are reported in l/min. When entering parameters for the gender specific equation for Z-scores, the parameter for both male and female must be converted to l/min. When entering parameters from reference papers the units used must be the same.

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